

## MCMBP

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**73kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human MCMBP

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

MCM-BP; C10orf119;

**Catalog #:**A3017

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**79892

**Isotype:**IgG

**Swiss Prot:**Q9BTE3

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a protein which is a component of the hexameric minichromosome maintenance (MCM) complex which regulates initiation and elongation of DNA. Multiple transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*