

GLRA3

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:200 - 1:1000 IHC 1:10 - 1:50

Calculated MW:54kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human GLRA3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A3024

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8001

Isotype:IgG

Swiss Prot:O75311

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the ligand-gated ion channel protein family. The encoded protein is a member of the glycine receptor subfamily. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).