

CDH1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IP

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IP 1:20 - 1:50

Calculated MW:135kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human CDH1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

dt

Synonym:

CDH1; Arc-1;CD324 ;CDHE ; ECAD;LCAM ; UVO ;E-Cad/CTF3 ;Uvomorulin; E-Cadherin;

Catalog #:A3044

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:999

Isotype:IgG

Swiss Prot:P12830

Purity:Affinity purification

For research use only.

Background:

This gene is a classical cadherin from the cadherin superfamily. The encoded protein is a calcium dependent cell-cell adhesion glycoprotein comprised of five extracellular cadherin repeats, a transmembrane region and a highly conserved cytoplasmic tail. Mutations in this gene are correlated with gastric, breast, colorectal, thyroid and ovarian cancer. Loss of function is thought to contribute to progression in cancer by increasing proliferation, invasion, and/or metastasis. The ectodomain of this protein mediates bacterial adhesion to mammalian cells and the cytoplasmic domain is required for internalization. Identified transcript variants arise from mutation at consensus splice sites. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).