

NCK2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:25 - 1:100

Calculated MW:43kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human NCK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GRB4; NCKbeta;

Background:

This gene encodes a member of the NCK family of adaptor proteins. The protein contains three SH3 domains and one SH2 domain. The protein has no known catalytic function but has been shown to bind and recruit various proteins involved in the regulation of receptor protein tyrosine kinases. It is through these regulatory activities that this protein is believed to be involved in cytoskeletal reorganization. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

To place an order, please [Click HERE](#).

Catalog #:A3047

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8440

Isotype:IgG

Swiss Prot:O43639

Purity:Affinity purification

For research use only.