

## AKAP8L

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:200 - 1:1000 IHC 1:25 - 1:100

**Calculated MW:**72kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human AKAP8L

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HA95; HAP95; NAKAP; NAKAP95;

*To place an order, please [Click HERE](#).*

**Catalog #:**A3050

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**26993

**Isotype:**IgG

**Swiss Prot:**Q9ULX6

**Purity:**Affinity purification

For research use only.