

## NeuN

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:200 - 1:1000 IHC 1:15 - 1:50

**Calculated MW:**34kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human NeuN

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FOX3; NEUN; FOX-3; HRNBP3;

*To place an order, please [Click HERE](#).*

**Catalog #:**A3051

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**146713

**Isotype:**IgG

**Swiss Prot:**A6NFN3

**Purity:**Affinity purification

For research use only.