

NCAN

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:25 - 1:100

Calculated MW:143kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human NCAN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CSPG3;

Catalog #:A3052

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1463

Isotype:IgG

Swiss Prot:O14594

Purity:Affinity purification

For research use only.

Background:

Neurocan is a chondroitin sulfate proteoglycan thought to be involved in the modulation of cell adhesion and migration.

To place an order, please [Click HERE](#).