

NTSR1

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:46kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human NTSR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NTR

Catalog #:A3054

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4923

Isotype:IgG

Swiss Prot:P30989

Purity:Affinity purification

For research use only.

Background:

Neurotensin receptor 1 belongs to the large superfamily of G-protein coupled receptors. NTSR1 mediates the multiple functions of neurotensin, such as hypotension, hyperglycemia, hypothermia, antinociception, and regulation of intestinal motility and secretion.

To place an order, please [Click HERE](#).