

## NTSR1

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**46kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human NTSR1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

NTR

**Catalog #:**A3054

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4923

**Isotype:**IgG

**Swiss Prot:**P30989

**Purity:**Affinity purification

For research use only.

**Background:**

Neurotensin receptor 1 belongs to the large superfamily of G-protein coupled receptors. NTSR1 mediates the multiple functions of neurotensin, such as hypotension, hyperglycemia, hypothermia, antinociception, and regulation of intestinal motility and secretion.

*To place an order, please [Click HERE](#).*