

GRIN2B

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:190kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human GRIN2B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

gi

Synonym:

MRD6; NR2B; hNR3; GluN2B; GRIN2B;

Catalog #:A3056

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2904

Isotype:IgG

Swiss Prot:Q13224

Purity:Affinity purification

For research use only.

Background:

N-methyl-D-aspartate (NMDA) receptors are a class of ionotropic glutamate receptors. NMDA receptor channel has been shown to be involved in long-term potentiation, an activity-dependent increase in the efficiency of synaptic transmission thought to underlie certain kinds of memory and learning. NMDA receptor channels are heteromers composed of three different subunits: NR1 (GRIN1), NR2 (GRIN2A, GRIN2B, GRIN2C, or GRIN2D) and NR3 (GRIN3A or GRIN3B). The NR2 subunit acts as the agonist binding site for glutamate. This receptor is the predominant excitatory neurotransmitter receptor in the mammalian brain.

To place an order, please [Click HERE](#).