

P2RY4

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:41kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human P2RY4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NRU; P2P; UNR; P2Y4

Catalog #:A3059

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5030

Isotype:IgG

Swiss Prot:P51582

Purity:Affinity purification

For research use only.

Background:

The product of this gene belongs to the family of G-protein coupled receptors. This family has several receptor subtypes with different pharmacological selectivity, which overlaps in some cases, for various adenosine and uridine nucleotides. This receptor is responsive to uridine nucleotides, partially responsive to ATP, and not responsive to ADP.

To place an order, please [Click HERE](#).