

PPAP2A

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:32kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human PPAP2A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LPP1; PAP2; LLP1a; PAP-2a;

Background:

The protein encoded by this gene is a member of the phosphatidic acid phosphatase (PAP) family. PAPs convert phosphatidic acid to diacylglycerol, and function in synthesis of glycerolipids and in phospholipase D-mediated signal transduction. This enzyme is an integral membrane glycoprotein that plays a role in the hydrolysis and uptake of lipids from extracellular space. Alternate splicing results in multiple transcript variants of this gene.

To place an order, please [Click HERE](#).

Catalog #:A3062

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8611

Isotype:IgG

Swiss Prot:O14494

Purity:Affinity purification

For research use only.