

PAX1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:55kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human PAX1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OFC2; HUP48;

Background:

This gene is a member of the paired box (PAX) family of transcription factors. Members of the PAX family typically contain a paired box domain and a paired-type homeodomain. These genes play critical roles during fetal development. This gene plays a role in pattern formation during embryogenesis and may be essential for development of the vertebral column. This gene is silenced by methylation in ovarian and cervical cancers and may be a tumor suppressor gene. Mutations in this gene are also associated with vertebral malformations.

To place an order, please [Click HERE](#).

Catalog #:A3066

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5075

Isotype:IgG

Swiss Prot:P15863

Purity:Affinity purification

For research use only.