

PAX2

Reactivity:Human Mouse Rat

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:200

Calculated MW:45kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human PAX2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

dk

Synonym:

Paired box 2

Catalog #:A3067

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5076

Isotype:IgG

Swiss Prot:Q02962

Purity:Affinity purification

For research use only.

Background:

This gene encodes type I adenylate cyclase activating polypeptide receptor, which is a membrane-associated protein and shares significant homology with members of the glucagon/secretin receptor family. This receptor mediates diverse biological actions of adenylate cyclase activating polypeptide 1 and is positively coupled to adenylate cyclase. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified.

To place an order, please [Click HERE](#).