

PIBF1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW:90kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human PIBF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

PIBF; CEP90; C13orf24; RP11-505F3.1

Catalog #:A3072

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10464

Isotype:IgG

Swiss Prot:Q8WXW3

Purity:Affinity purification

For research use only.

Background:

Progesterone-induced-blocking factor 1 is a protein that in humans is encoded by the PIBF1 gene. It has been shown to localize to the centrosome and has also been named CEP90. Mediator of progesterone that by acting on the phospholipase A2 enzyme interferes with arachidonic acid metabolism, induces a Th2 biased immune response, and by controlling NK activity exerts an anti-abortive effect.

To place an order, please [Click HERE](#).