

## PFN2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**15kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human PFN2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

PFL; D3S1319E;

**Background:**

The protein encoded by this gene is a ubiquitous actin monomer-binding protein belonging to the profilin family. It is thought to regulate actin polymerization in response to extracellular signals.

There are two alternatively spliced transcript variants encoding different isoforms described for this gene.

*To place an order, please [Click HERE](#).*

**Catalog #:**A3074

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5217

**Isotype:**IgG

**Swiss Prot:**P35080

**Purity:**Affinity purification

For research use only.