

RNF126

Reactivity:Human Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:25 - 1:100

Calculated MW:36kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human RNF126

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A3081

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:55658

Isotype:IgG

Swiss Prot:Q9BV68

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a RING finger domain, a motif present in a variety of functionally distinct proteins and known to be involved in protein-protein and protein-DNA interactions.

To place an order, please [Click HERE](#).