

TCF7

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:42kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human TCF7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

kt

Synonym:

TCF-1

Catalog #:A3091

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6932

Isotype:IgG

Swiss Prot:P36402

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a transcriptional activator that plays an important role in lymphocyte differentiation. This gene is expressed predominantly in T-cells. The encoded protein can bind an enhancer element and activate the CD3E gene, and it also may repress the CTNNB1 and TCF7L2 genes through a feedback mechanism. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).