

TRAF7

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW:67kD

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human TRAF7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1ml

Synonym:

RFWD1; RNF119;

Catalog #:A3095

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:84231

Isotype:IgG

Swiss Prot:Q6Q0C0

Purity:Affinity purification

For research use only.

Background:

Tumor necrosis factor (TNF; see MIM 191160) receptor-associated factors, such as TRAF7, are signal transducers for members of the TNF receptor superfamily (see MIM 191190). TRAFs are composed of an N-terminal cysteine/histidine-rich region containing zinc RING and/or zinc finger motifs; a coiled-coil (leucine zipper) motif; and a homologous region that defines the TRAF family, the TRAF domain, which is involved in self-association and receptor binding.

To place an order, please [Click HERE](#).