

NEUROG1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW:26kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human NEUROG1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

AKA; ngn1; Math4C; bHLHa6; NEUROD3

Background:

Neurogenin-1

To place an order, please [Click HERE](#).

Catalog #:A3110

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4762

Isotype:IgG

Swiss Prot:Q92886

Purity:Affinity purification

For research use only.