

## ENO2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**47kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human ENO2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

NSE

**Catalog #:**A3118

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2026

**Isotype:**IgG

**Swiss Prot:**P09104

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes one of the three enolase isoenzymes found in mammals. This isoenzyme, a homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development in rats and primates.

*To place an order, please [Click HERE](#).*