

## PPP2CB

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**36kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human PPP2CB

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

qt

**Synonym:**

PP2CB; PP2Abeta;

**Catalog #:**A3122

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5516

**Isotype:**IgG

**Swiss Prot:**P62714

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes the phosphatase 2A catalytic subunit. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. This gene encodes a beta isoform of the catalytic subunit.

**To place an order, please [Click HERE](#).**