

## PTN

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:20 - 1:100

**Calculated MW:**19kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic Peptide of human PTN

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HARP; HBNF; HBGF8; NEGF1

**Background:**

Pleiotrophin

*To place an order, please [Click HERE](#).*

**Catalog #:**A3124

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5764

**Isotype:**IgG

**Swiss Prot:**P21246

**Purity:**Affinity purification

For research use only.