

SCTR

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:50kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SCTR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SR

Catalog #:A3125

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6344

Isotype:IgG

Swiss Prot:P47872

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a G protein-coupled receptor and belongs to the glucagon-VIP-secretin receptor family. It binds secretin which is the most potent regulator of pancreatic bicarbonate, electrolyte and volume secretion. Secretin and its receptor are suggested to be involved in pancreatic cancer and autism.

To place an order, please [Click HERE](#).