

SLC2A4RG

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:41kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human SLC2A4RG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GEF; HDBP1; HDBP-1; Si-1-2; Si-1-2-19;

Background:

The protein encoded by this gene is a nuclear transcription factor involved in the activation of the solute carrier family 2 member 4 gene. The encoded protein interacts with another transcription factor, myocyte enhancer factor 2, to activate transcription of this gene.

To place an order, please [Click HERE](#).

Catalog #:A3127

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:56731

Isotype:IgG

Swiss Prot:Q9NR83

Purity:Affinity purification

For research use only.