

## SLC22A17

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:25 - 1:100

**Calculated MW:**58kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human SLC22A17

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

BOCT; BOIT; 24p3R; NGALR; hBOIT; NGALR2; NGALR3;

*To place an order, please [Click HERE](#).*

**Catalog #:**A3128

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**51310

**Isotype:**IgG

**Swiss Prot:**Q8WUG5

**Purity:**Affinity purification

For research use only.