

## MLF1

---

**Reactivity:**Human Mouse

**Tested applications:**WB IHC IP

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

**Calculated MW:**31kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human MLF1

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:**A3960

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4291

**Isotype:**IgG

**Swiss Prot:**P58340

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes an oncoprotein which is thought to play a role in the phenotypic determination of hemopoetic cells. Translocations between this gene and nucleophosmin have been associated with myelodysplastic syndrome and acute myeloid leukemia. Multiple transcript variants encoding different isoforms have been found for this gene.

**To place an order, please [Click HERE](#).**