

## MOV10

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF IP

**Recommended Dilution:**WB 1:1000 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:100 IP 1:20 - 1:50

**Calculated MW:**114kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human MOV10

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

gb110; fSAP113;

*To place an order, please [Click HERE](#).*

**Catalog #:**A3966

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4343

**Isotype:**IgG

**Swiss Prot:**Q9HCE1

**Purity:**Affinity purification

For research use only.