

NDUFB9

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:22kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human NDUFB9

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

i

Synonym:

B22; LYRM3; CI-B22; UQOR22;

To place an order, please [Click HERE](#).

Catalog #:A3981

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4715

Isotype:IgG

Swiss Prot:Q9Y6M9

Purity:Affinity purification

For research use only.