

## NFE2L1

---

**Reactivity:**Human Mouse

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:200 - 1:2000 IHC 1:50- 1:200

**Calculated MW:**85kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human NFE2L1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:**A3988

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4779

**Isotype:**IgG

**Swiss Prot:**Q14494

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a protein that is involved in globin gene expression in erythrocytes. Confusion has occurred in bibliographic databases due to the shared symbol of NRF1 for this gene, NFE2L1, and for "nuclear respiratory factor 1" which has an official symbol of NRF1. [provided by RefSeq, Jul 2008]

*To place an order, please [Click HERE](#).*