

NVL

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:73kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NVL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Catalog #:A3995

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4931

Isotype:IgG

Swiss Prot:O15381

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the AAA (ATPases associated with diverse cellular activities) superfamily. Multiple transcript variants encoding different isoforms have been found for this gene. Two encoded proteins, described as major and minor isoforms, have been localized to distinct regions of the nucleus. The largest encoded protein (major isoform) has been localized to the nucleolus and shown to participate in ribosome biosynthesis (PMID: 15469983, 16782053), while the minor isoform has been localized to the nucleoplasmin.

To place an order, please [Click HERE](#).