

## OSBP

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 89kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human OSBP

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

OSBP1;

**Catalog #:** A3998

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 5007

**Isotype:** IgG

**Swiss Prot:** P22059

**Purity:** Affinity purification

For research use only.

**Background:**

Oxysterol binding protein is an intracellular protein that is believed to transport sterols from lysosomes to the nucleus where the sterol down-regulates the genes for the LDL receptor, HMG-CoA reductase, and HMG synthetase

*To place an order, please [Click HERE](#).*