

## PAFAH1B2

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**26kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human PAFAH1B2

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HEL-S-303

**Catalog #:**A4001

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5049

**Isotype:**IgG

**Swiss Prot:**P68402

**Purity:**Affinity purification

For research use only.

**Background:**

Platelet-activating factor acetylhydrolase (PAFAH) inactivates platelet-activating factor (PAF) into acetate and LYSO-PAF. This gene encodes the beta subunit of PAFAH, the other subunits are alpha and gamma. Multiple alternatively spliced transcript variants have been described for this gene.

*To place an order, please [Click HERE](#).*