

## PGAM1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:2000

**Calculated MW:**29kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human PGAM1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

PGAMA; PGAM-B;

**Catalog #:**A4015

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5223

**Isotype:**IgG

**Swiss Prot:**P18669

**Purity:**Affinity purification

For research use only.

**Background:**

Phosphoglyceric acid mutase (EC 2.7.5.3) is widely distributed in mammalian tissues where it catalyzes the reversible reaction of 3-phosphoglycerate (3-PGA) to 2-phosphoglycerate (2-PGA) in the glycolytic pathway (summary by Chen et al., 1974 [PubMed 4811757]).

*To place an order, please [Click HERE](#).*