

## SERPINA4

---

**Reactivity:**Human

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:20 - 1:200

**Calculated MW:**49kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human SERPINA4

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

KAL; KST; PI4; KLST; PI-4; kallistatin;

**To place an order, please [Click HERE](#).**

**Catalog #:**A4019

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5267

**Isotype:**IgG

**Swiss Prot:**P29622

**Purity:**Affinity purification

For research use only.