

## PIK3C3

---

**Reactivity:**Human

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:200 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:**100kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human PIK3C3

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

cq

**Synonym:**

VPS34; hVps34;

**To place an order, please [Click HERE](#).**

**Catalog #:**A4021

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5289

**Isotype:**IgG

**Swiss Prot:**Q8NEB9

**Purity:**Affinity purification

For research use only.