

PPP6C

Reactivity: Human Mouse Rat

Tested applications: WB IP

Recommended Dilution: WB 1:200 - 1:2000 IP 1:20 - 1:50

Calculated MW: 35kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PPP6C

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

e

Synonym:

PP6; PP6C;

Catalog #: A4039

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5537

Isotype: IgG

Swiss Prot: O00743

Purity: Affinity purification

For research use only.

Background:

This gene encodes the catalytic subunit of protein phosphatase, a component of a signaling pathway regulating cell cycle progression. Splice variants encoding different protein isoforms exist. The pseudogene of this gene is located on chromosome X.

To place an order, please [Click HERE](#).