

## PSMD5

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IP

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

**Calculated MW:**56kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human PSMD5

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

S5B;

**Catalog #:**A4056

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5711

**Isotype:**IgG

**Swiss Prot:**Q16401

**Purity:**Affinity purification

For research use only.

**Background:**

The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure composed of 2 complexes, a 20S core and a 19S regulator. The 20S core is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. The 19S regulator is composed of a base, which contains 6 ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-ATPase subunits.

Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. This gene encodes a non-ATPase subunit of the 19S regulator base that functions as a chaperone protein during 26S proteasome assembly.

**To place an order, please [Click HERE](#).**