

## RBBP5

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:200 - 1:2000

**Calculated MW:**59kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human RBBP5

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

RBQ3; SWD1;

**Catalog #:**A4072

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5929

**Isotype:**IgG

**Swiss Prot:**Q15291

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a ubiquitously expressed nuclear protein which belongs to a highly conserved subfamily of WD-repeat proteins. The encoded protein binds directly to retinoblastoma protein, which regulates cell proliferation. It interacts preferentially with the underphosphorylated retinoblastoma protein via the E1A-binding pocket B. Three alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.

*To place an order, please [Click HERE](#).*