

RBM4

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:40kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic Peptide of human RBM4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

LARK; RBM4A; ZCRB3A; ZCCHC21;

To place an order, please [Click HERE](#).

Catalog #:A4073

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5936

Isotype:IgG

Swiss Prot:Q9BWF3

Purity:Affinity purification

For research use only.