

RNH1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 50kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RNH1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

a

Synonym:

RAI; RNH;

Catalog #: A4079

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6050

Isotype: IgG

Swiss Prot: P13489

Purity: Affinity purification

For research use only.

Background:

Placental ribonuclease inhibitor (PRI) is a member of a family of proteinaceous cytoplasmic RNase inhibitors that occur in many tissues and bind to both intracellular and extracellular RNases (summarized by Lee et al., 1988)

To place an order, please [Click HERE](#).