

RNPEP

Reactivity:Human Mouse Rat

Tested applications:WB IF

Recommended Dilution:WB 1:1000 - 1:2000 IF 1:10 - 1:100

Calculated MW:73kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RNPEP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A4080

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6051

Isotype:IgG

Swiss Prot:Q9H4A4

Purity:Affinity purification

For research use only.