

RTN2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW:59kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human RTN2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NSP2; NSPL1; NSPLI; SPG12;

Catalog #:A4100

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6253

Isotype:IgG

Swiss Prot:O75298

Purity:Affinity purification

For research use only.

Background:

This gene belongs to the family of reticulon encoding genes. Reticulons are necessary for proper generation of tubular endoplasmic reticulum and likely play a role in intracellular vesicular transport. Alternatively spliced transcript variants encoding different isoforms have been identified. Mutations at this locus have been associated with autosomal dominant spastic paraplegia-12.

To place an order, please [Click HERE](#).