

S100A6

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 10kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human S100A6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

2A9; PRA; 5B10; CABP; CACY;

Catalog #: A4104

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6277

Isotype: IgG

Swiss Prot: P84098

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein may function in stimulation of Ca²⁺-dependent insulin release, stimulation of prolactin secretion, and exocytosis. Chromosomal rearrangements and altered expression of this gene have been implicated in melanoma.

To place an order, please [Click HERE](#).