

SH3GL3

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:4000

Calculated MW:39kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SH3GL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CNSA3; EEN-B2; SH3D2C; SH3P13; HsT19371;

To place an order, please [Click HERE](#).

Catalog #:A4112

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6457

Isotype:IgG

Swiss Prot:Q99963

Purity:Affinity purification

For research use only.