

## SNRPG

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:200 - 1:2000

**Calculated MW:**8kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human SNRPG

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

SMG; Sm-G;

**To place an order, please [Click HERE](#).**

**Catalog #:**A4122

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**6637

**Isotype:**IgG

**Swiss Prot:**P62308

**Purity:**Affinity purification

For research use only.