

TROVE2

Reactivity:Human Mouse Rat

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:10 - 1:100

Calculated MW:61kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human TROVE2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RO60; SSA2; RORNP;

To place an order, please [Click HERE](#).

Catalog #:A4128

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6738

Isotype:IgG

Swiss Prot:P10155

Purity:Affinity purification

For research use only.