

## TBX6

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IP

**Recommended Dilution:** WB 1:500 - 1:2000 IP 1:20 - 1:50

**Calculated MW:** 47kDa

**Observed MW:** Refer to figures

**Immunogen:**

A synthetic Peptide of human TBX6

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

k

**Synonym:**

SCDO5;

**Catalog #:** A4142

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 6911

**Isotype:** IgG

**Swiss Prot:** O95947

**Purity:** Affinity purification

For research use only.

**Background:**

This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental processes. Knockout studies in mice indicate that this gene is important for specification of paraxial mesoderm structures.

*To place an order, please [Click HERE](#).*