

## TEAD4

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:200 - 1:2000

**Calculated MW:**48kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human TEAD4

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

kn

**Synonym:**

TEF3; RTEF1; TEF-3; EFTR-2; TEFR-1; TCF13L1; hRTEF-1B;

**Catalog #:**A4151

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**7004

**Isotype:**IgG

**Swiss Prot:**Q15561

**Purity:**Affinity purification

For research use only.

**Background:**

This gene product is a member of the transcriptional enhancer factor (TEF) family of transcription factors, which contain the TEA/ATTS DNA-binding domain. It is preferentially expressed in the skeletal muscle, and binds to the M-CAT regulatory element found in promoters of muscle-specific genes to direct their gene expression. Alternatively spliced transcripts encoding distinct isoforms, some of which are translated through the use of a non-AUG (UUG) initiation codon, have been described for this gene.

**To place an order, please [Click HERE](#).**