

## UBP1

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:2000

**Calculated MW:**60kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human UBP1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

k

**Synonym:**

LBP1A; LBP1B; LBP-1B; LBP-1a;

*To place an order, please [Click HERE](#).*

**Catalog #:**A4177

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**7342

**Isotype:**IgG

**Swiss Prot:**Q9NZI7

**Purity:**Affinity purification

For research use only.